

## PRO SPIKE

Compact container for used needles



Disposable container for collection and dispose of used needles.

### Specific features

Small and rigid structure

It has a clever lid with a safe closing system and a clip so you can attach it to belts, backpacks and other things

The kit includes 5 containers

### Technical data

Dimensions	90x54x185 mm
Weight <sup>(one piece)</sup>	88 g

PRO SPIKE BOX 5 PCS - COMPACT  
CONTAINER FOR USED NEEDLES  
**FG04591C**

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)