

Collection canister



Collection canister for small volumes intended for “High vacuum – High flow” applications. The device can be used with vacuum regulators or Venturi suction devices.

Specific Characteristics

Autoclavable canister
Overfill floating valve to ensure maximum safety against contamination

Technical data

Volume	300 ml
Dimensions	180x98x70 mm
Supply vacuum	Max -950mbar / 5min
Overfill valve	Floating valve
Inlet connection	ISO G 1/2" F
Outlet connection	Hose \varnothing 8,0 ÷ 9,2 mm
Autoclavable	121°C 15 min
Weight	180 g

Class I Med.Device according to Dir. 93/42/CEE

300cc VASE WITH CAP FOR
VENTURI ASPIRATOR

FL08221E

CND Classification A0699

Registration number 154294



Distributed by Spencer Italia s.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it