

FL-X Humidifier



Humidifier intended to be connected to a flowmeter within an oxygen distribution systems.

Specific features

The humidifier helps to reduce the dryness of mucosas during administration of oxygen

Hose output for the connection of devices used for oxygen administration

Technical data

Dimensions	ø 70x135 mm
Max. volume	130 ml
Output hose	ø 7 mm
Operating pressure	350 kPa
Max pressure	500 kPa
Max flow	20 L/min
Materials	PC, Al
Weight	110 g

Compliant with

UNI EN ISO 8185

Class IIa Medical device according to dir. 93/42/CE

FL-X HUMIDIFIER ONLY
FL08302B



Registration number 21663

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it