

## FL-X

### Flowmeter with UNI connection probe



Flowmeter with a male probe intended to be connected to gas-specific sockets compliant with UNI standards. The socket shall be installed within a distribution system for medical gases.

#### Specific Characteristics

Visualization of the gas flow delivered to the patient  
 Safety and rapidity of the connection and disconnection procedures, thanks to the gas-specific probe  
 Hose holder output allow quick connection of oxygen tubes

#### Technical data

Dimensions	133x177x40 mm
Nominal pressure	350 kPa
Flow	From 0 to 15 l/min
Deliverable gas	Oxygen
Materials	PC, anodized Al, Brass
Weight	245 g
Output	Hose holder $\varnothing$ 7mm

#### Compliant with

UNI EN 15002  
 UNI 9507

Medical device CE marked according to dir. 93/42/CE

FL-X FLOW. W/UNI  
 CONNECTOR AND HOSE  
 HOLDER  $\varnothing$  7

**FL08320C**

**Registration number** 24304

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)