

FL-X

Flowmeter with DIN connection probe



Flowmeter with a male probe intended to be connected to gas-specific sockets compliant with DIN standards. The socket shall be installed within a distribution system for medical gases.

Specific Characteristics

Visualization of the gas flow delivered to the patient
 Safety and rapidity of the connection and disconnection procedures, thanks to the gas-specific probe
 Hose holder output allow quick connection of oxygen tubes

Technical data

Dimensions	120x177x40 mm
Nominal pressure	350 kPa
Flow	From 0 to 15 l/min
Deliverable gas	Oxygen
Materials	PC, anodized Al, Brass
Weight	200 g
Output	Hose holder \varnothing 7mm

Compliant with

UNI EN 15002
 DIN 13260-2

Medical device CE marked according to dir. 93/42/CE

FL-X FLOW W/DIN CONNECTOR
 AND HOSE HOLDER \varnothing 7

FL08322C

Registration number 24343

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it