

MATRIX THJ

Holster



Holster characterized by compact dimensions and compartments designed to store laryngoscope or torch.

Specific features

- Semi-rigid durable structure
- Strap to attach the holster to the belt
- Elastic bands to store accessories
- Two mesh pockets on the internal side of the two halves
- Closure with double zip

Technical data¹

Dimensions	210x45x60 ± 10mm
Weight	200 ± 10 g
Material	PU

HOLSTER MATRIX ORANGE
GD01002E

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it