

## RING CUTTER



Device for manual cutting of rings.

### Specific features

Sterilizable  
Toothed blade to increase effectiveness of the cut  
The size of the handle and its design, increase stability during cutting

### Technical data

Length	165 mm
Blade diameter	20,5 mm
Weight	91 g
Handle length	100 mm
Structure material	Chromed brass
Blade material	Steel
Cutting direction	Clockwise
Cutting zone width	17 mm

Class I Med.Device according to Dir. 93/42/CEE

RING CUTTER - STEEL RING  
CUTTER W/BRASS HANDLE  
GD02020E



Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l a socio unico. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)