

Match 1

Glass breaking pen



Device to brake vehicle's glasses during on-road rescue operations .

Specific Characteristics

Extremely compact
Knurled surface

Technical data

Dimensions	ø12x124 mm
Material	Steel
Weight	60 g

MATCH 1 - GLASS BREAKING
PEN
GD02030E

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it