

# MATCH 2 COMPACT

## Emergency hammer



Emergency hammer used for window breaking and seat belt cutting.

### Specific features

The hammer is supplied with a wall fastener

Easy and intuitive to use

### Technical data

Hammer dimensions	135x75x25 mm
Dimensions with fastener	140x90x35 mm
Hammer weight	150g
Overall weight	160 g
Material	ABS, stainless steel

### Spare parts

GD02041 Wall fastener

Class I Med.Device according to Dir. 93/42/CEE

MATCH 2 COMPACT – EMERGENCY  
HAMMER

**GD02040C**

Rev.1 (01/04/2019)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)