

MATCH 3

Multi-purpose scissors



1 Used to break glass.



2 Used to cut seatbelts.



3 Used to cut clothes.

Multi-purpose scissors with wall mount to provide safe installation inside any vehicle.

Specific features

- Interlocking fastening allows a quick access
- The support ensures that the device remains closed when not used
- The fastener is installed by means of 3 screws for maximum safety
- Rounded tip to limit risks of injury to the skin
- The shape of the plastic part covering the blade, allows to grasp the scissors using them as an evacuation hammer

Technical data

Scissors dimensions	247x105x18 mm
Dimensions with fastener	252x115x22 mm
Blade length	100 mm
Material	Steel, Nylon
Scissors weigh	374 g
Fastener weight	36 g

Spare parts

- GD02064 Wall fastener
- GD02069 Hooks for carpet fastening

Class I Med.Device according to Dir. 93/42/CEE

MATCH 3 MULTIFUNCTIONAL
SCISSORS

GD02060C

CND Classification L01040199

Registration number 161376



NATO STOCK N° 6515-15-010-7023

Rev.3 (07/06/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it