

## Matrix MHP

Holster for multi-purpose scissors



Holster made of thermoformed material deal to store the multi-purpose scissors Spencer Match 3.

### Specific features

The semi-rigid structure limits risks for the operator and provides excellent protection to the content

Internal mesh pockets for accessories such as gloves or gauzes

Strap for belt attachment

Closure with double zip for easy access

### Technical data

Dimensions	255x125x45 mm
Materials	EVA, PU
Weight	120 g

MATRIX MHP YELLOW POLYURETHANE – HOLSTER  
FOR MULTI-PURPOSE SCISSORS  
GD02071E

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)