

Matrix MHJ

Holster for multi-purpose scissors



Holster made of thermoformed material deal to store the multi-purpose scissors Spencer Match 3.

Specific features

The semi-rigid structure limits risks for the operator and provides excellent protection to the content

Internal mesh pockets for accessories such as gloves or gauzes

Strap for belt attachment

Closure with double zip for easy access

Technical data

Dimensions	255x125x45 mm
Materials	EVA, Jersey
Weight	120 g

MATRIX MHJ ORANGE

GD02072E

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it