

## Eye medication KIT



Eye medication kit intended for first aid cleaning and bandage.

Content of the kit	
Article	Q.ty
Eye washing solution	1
Hypoallergenic non-sterile patch	1
Non-stick sterile pad with adhesive edge	3
Sterile eye pad with adhesive edge	2
First aid short guide	1

Technical data	
Eye solution	100 ml NaCl 0,9% Sterile
Non-stick pad	50x70 mm
Eye pad	60x45 mm circa
Patch	2,5 cm x 5m
Weight of the kit	220 g

MD CE marked according to dir. 93/42/CE

EYE MEDICATION KIT

**GM00109K**

Distributed by Spencer Italia srl

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)