

GMX 234

Curved scissors with rounded end



Stainless steel curved scissors with rounded end.

Specific features

Sterilizable
Satin finish to reduce light reflection

Technical data

| | |
|--------------|-----------------|
| Total length | 161 mm |
| Blade length | 45 mm |
| Weight | 55 g |
| Material | Stainless steel |

Class I Med.Device according to Dir. 93/42/CEE

SURGICAL CURVED SCISSORS
ROUND END 16cm

GM00234G

CND Classification



L031306

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it