

## GMX 254

### Curved scissors



Stainless steel curved scissors with pointed end.

#### Specific features

Sterilizable  
Satin finish to reduce light reflection

#### Technical data

Total length	160 mm
Blade length	45 mm
Weight	52 g
Material	Stainless steel

Class I Med.Device according to Dir. 93/42/CEE

SURGICAL CURVED SCISSORS  
POINTED END 16 cm

**GM02254G**

**CND Classification**



L031306

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l a socio unico. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)