

PLASTIC KIDNEY DISH



Plastic kidney dish.

Specific features

Polypropylene body
Glossy surface

Technical data

Dimensions	190x79x38 mm
Material	PP
Weight	48 g

Class I Med. Dev. According to Dir. 93/42/CEE

KIDNEY DISH PLASTIC
GM00283G



Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it