

GMX 363

Hemostatic forceps



Klemmer forceps, commonly used as hemostat for small blood vessels.

Specific features

Sterilizable
 Reeded distal part to offer a better grip.
 The absence of screws, makes it very stable and reliable

Technical data

Total length	158 mm
Clamping and reeded area length	45 mm
Weight	40 g
Material	Stainless steel

Class I Med.Device according to Dir. 93/42/CEE

GMX 363 - KLEMMER
 FORCEPS

GM00363A



CND Classification L031306

Registration number 135896

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it