

GMX 373

Hemostatic forceps



Klemmer curved forceps, commonly used as hemostat for small blood vessels.

Specific features

- Sterilizable
- Reeded distal part to offer a better grip.
- The absence of screws, makes it very stable and reliable

Technical data

Total length	158 mm
Clamping and reeded area length	45 mm
Weight	34 g
Material	Stainless steel

Class I Med.Device according to Dir. 93/42/CEE

GMX 373 - KLEMMER
HEAMOSTATIC CURVED
TWEEZER



GM00373A

CND Classification	L031306
Registration number	135908

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it