

GMX 412

Hemostatic forceps



Halsted-Mosquito forceps, commonly used as hemostat for small blood vessels.

Specific features

- Sterilizable
- Reeded distal part to offer a better grip.
- The absence of screws, makes it very stable and reliable

Technical data

Total length	125 mm
Clamping and reeded area length	22 mm
Weight	19 g
Material	Stainless steel

Class I Med.Device according to Dir. 93/42/CEE

GMX 412 - HAEMOST.
FORCEPS HALSTED-
MOSQUITO 123 mm

GM00412A



CND Classification	L031306
Registration number	135914

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l a socio unico. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it