

GMX 563

Mayo-Hegar needle holder



Mayo-Hegar needle driver to hold a suturing needle during suturing and surgical procedures.

Specific features

The surface finish doesn't reflect light during use
Particularly useful in deep sutures where there is not enough space for fingers

Technical data

Length	14 cm
Material	Steel
Weight	32 g

Class I MD according to Dir. 93/42/CEE

GMX 563 - MAYO HEGAR NEEDLE
HOLDER 14CM
GM00563A



CND Classification	L0205
Registration number	135996

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it