

S-G

Sterile gauze



Sterile hydrophilic gauze for medication.

Specific features

25 sterile disposable gauzes
Sealed envelope

Technical data

Single gauze dimensions	10x10 cm
Pachage weight	7 g
Gauze material	100% cotton
Sterilization method	EO

Class I S Med. Dev. According to Dir. 93/42/CEE

GAUZE COMPRESSES STERILE
10X10CM 25 SHEETS PER PAC
GM02101K



Registration number	180163/R
CND Classification	M0201010101

Distributed by Spencer Italia S.r.l

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it