

S-G Sterile Gauze



Sterile gauze for dressings.

Specific Characteristics

Sealed package
Disposable

Technical data

Dimensions opened gauze 18 x 40 cm

Class Is MD according to Dir. 93/42/CEE

STERIL GAUZE PAD 18X40CM

GM02103K



Distributed by Spencer Italia S.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it