

# T-plaster

## Roll of plaster



Adhesive patch with zinc oxide.

### Specific Characteristics

Particularly suitable to the treatment of the lower limbs  
 Hypoallergenic  
 Soothing action  
 Can be cut without the use of scissors

### Technical data

Dimensions	5cm x 5 m
Weight	60 g

**Class I MD according to Dir. 93/42/CEE**

CLOTH PLASTER 5CM X 5M  
 GM2250K



Distributed by Spencer Italia S.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)