

## NET-B

### Net bandage



Elastic tubular bandage suitable for bandages and dressings. The bandage should not be applied on broken skin.

#### Specific features

Estremamente elastica  
Possibility of sterilization  
Particularly suitable for application on fingers

#### Technical data

Stretched length	3 m
Caliber	0
Weight (bandage and box)	15 g
Bandage material	27% Latex, 73% polyamide
Sterilizable by	Autoclave 121°, Gamma ray

#### Class I MD according to Dir. 93/42/CEE

NET-B - BANDAGE FOR FINGER  
GM03000K



CND Classification

M03030103

Distributed by Spencer Italia S.r.l

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)