

## GMX 02

Sterile disposable syringe 2,5ml



Disposable sterile disposable syringe.

### Specific features

EO Sterilized

Shelf life 5 year with intact package

### Technical data

Useful capacity	2,75 ml
Scale division	0,1 ml
Numbers on scale	Every 0,5 ml
Tolerance on the supplied capacity	1,5ml $\pm$ 5%
Seal capacity (compression)	Min 300 kPa
Death space	0,07 ml

Class IIa med.dev. according to dir. 93/42/CEE

STERILE DISPOSABLE SYRINGE

2,5ML

GM03020K



Rev.0 (28/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)