

GMX 05

Sterile disposable syringe 5ml



Disposable sterile disposable syringe.

Specific features

EO Sterilized
Shelf life 5 year with intact package

Technical data

Useful capacity	5,5 ml
Scale division	0,2 ml
Numbers on scale	Every 1 ml
Tolerance on the supplied capacity	2,5ml \pm 4%
Seal capacity (compression)	Min 300 kPa
Death space	0,07 ml

Class IIa med.dev. according to dir. 93/42/CEE

STERILE DISPOSABLE SYRINGE 5ML

GM03050K



Rev.0 (28/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it