

S-plaster

Roll of plaster



Not sterile plaster roll, to be used for dressing in non-invasive application on intact skin.

Specific Characteristics

Indicated to secure any type of dressing such as gauzes, probes, catheters and electrodes

The hypoallergenic artificial silk improves transpiration, allows easy application and makes the plaster imperceptible to the patient

Delicate upon removal, doesn't leave adhesive residues on the skin

Can be easily cut without the use of scissors

Technical data

Dimensions	1,25 cm x 9,14 m
Weight	16 g

Class I MD according to Dir. 93/42/CEE

S-PLASTER - SILK BANDAID 1,25cm x 9,14m

GM03215K

CND Classification



M050101

Distributed by Spencer Italia S.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it