

SSP Plasters



Plasters to keep wounds in best conditions for rapid healing.

Specific features

The perforated support makes the patch soft and elastic
Single package

Technical data

The box contains:	5 strips 7x2 cm
	10 strips 6x2 cm
	5 strips 4x1 cm
Box weight	15 g

Class I Med.Device according to Dir. 93/42/CEE

STRIP BANDAGES 20 PCS
DIFFERENT SIZES

GM04021K

CND Classification	M040199
Registration number	207649

Distributed by Spencer Italia s.r.l

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it