

SPENCER SICK SACK

Vomit bag



12 Transparent vomit bag.

Specific features

- Transparent to allow content view
- Carton inlet to provide rigidity
- Graduated bag
- Membrane to prevent spillage

Technical data

Package dimensions (12 pcs)	135x194x20 mm
Package dimensions (1 piece)	135x194 mm
Bag dimensions	390x180
Bag volume	1000cc
Weight (1 piece)	20 g
Weight of the kit	230 g

12 VOMIT BAGS

GM22222C

Rev.1 (25/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it