

## TIV Automatic tourniquet





Tourniquet equipped with a mechanism allowing gradual or instant loosening.

**Specific Characteristics** 

- Mechanism is simple to operate
- The elastic is large to act on a wider surface

Technical data	
Elastic width	26 mm
Elastic length <sup>(not extended)</sup>	From 0 to 460 mm
Weight	32 g

Class I Med. Device according to Dir. 93/42/CEE	
T.I.V. SPENCER AUTOMATIC	
TOURNIQUETTE	
IF00100C	
CND Classification	M040599
Registration number	270269

Rev.2 (15/01/2018) UNCHECKED COPY – further revisions will be available on http://support.spencer.it Spencer Italia S.r.I. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222