

# IFX50

## Tourniquet



Tubular Tourniquet.

### Specific features

Traditional type Tubular tourniquet

### Technical data

Length	38 cm
Weight	6 g
Material	Latex rubber

Class I Med.Device according to Dir. 93/42/CEE

IFX50 - RUBBER TOURNIQUETTE

IF00150C



CND Classification	V9003
Registration number	183276

Distributed by Spencer Italia s.r.l.

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)