

WIV

IV subministration bag



The WIV infusion subministration bag is a device that allows pressure to be applied to the infusion bag in order to empty the bag inside.

Specific features

- The transparent pocket used to store IV infusion allows a quick view of the content
- Eyelet on the upper part allow to hang the bag
- The integrated pressure gauge and the control valve allow extremely accurate pressure settings
- The device can be used with 500 and 1000 ml IV bags

Technical data

Dimensions	370 x 180 ± 10 mm
Length of pressure gauge tube	190 ± 10 mm
Length of inflation tube	190 ± 10 mm
Pressure gauge scale	0-300 mmHg Division every 2 mmHg
Maximum applicable pressure	250 mmHg
Pressure gauge diameter	50 mm
Materials	Nylon, PVC, Al, Chromed-plated brass
Weight	320 ± 20 g

Class I MD compliant with UE Reg. 2017/745

PRESSURE INFUSER WIV

IF00301A



CND Classification	V9099
Registration number	103738

Rev.0 (23/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it