

WIV

IV subministration Bag



Reusable device intended for IV infusion treatment, to be used during transportation of patient inside healthcare environments like hospitals, clinics and medical surgeries.

Specific features

- The transparent pocket used to store IV infusion, allows a quick view of the content
- Eyelet on the upper part, allow to hang the bag
- The integrated pressure gauge and the control valve, allow extremely accurate pressure settings
- The device can be used with 500cc and 1000cc IV bags

Technical data

Dimensions	370x180 ± 10mm
Press.gauge tube length	190 ± 10mm
Inflation bulb tube length	190 ± 10mm
Pressure gauge scale	0-300 mmHg Divisions every 2 mmhg
Maximum applicable pressure	250 mmHg
Pressure gauge diameter	50 mm
Materials	Nylon, rubber, Al, Chromed brass
Weight	320 ± 20g

Class I Med.Device according to Dir. 93/42/CEE

PRESSURE INFUSER WIV
IF00301A



CND Classification	V9099
Registration number	103738

Rev.3 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it