

Ice Pack

Instanti ce in PE bag



Instant ice pack to provide relief from pain in the event of trauma that requires cryotherapy and the use of cold.

Specific features

PE bag
 Quick cooling up to -5°C ($\pm 3^{\circ}\text{C}$)
 Duration of the cooling action approx 30 min
 25 disposable bags

Technical data

Dimensions	135x170 mm
Weight (1 piece)	230 g
Substance contained	Nitrogen salts
Bag material	PE

Class IIa med. Dev. according to Dir. 93/42/CEE

INSTANT ICE PE BAG 25PCS

IF00430A



CND Classification

M9001

Distributed from Spencer Italia s.r.l.

Rev.1 (15/01/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it