

## TRACK 4

### Telescopic IV pole



Track 4 is a telescopic IV pole for Spencer stretcher, used to maintain the hydraulic head needed to carry out infusion also during transport of the patient on the stretcher.

#### Specific features

- Easy to install
- Robust and lightweight structure
- It can be folded to minimize dimensions
- The locking rings allow a quick and effective adjustment
- Can hold two infusion bags
- Bands with straps to prevent oscillation of the bags during transport

#### Compatible stretchers

- Stretcher 506
- Cross Chair – Stretcher/Chair
- Stretcher Cross Up 8409

#### Technical data <sup>1</sup>

Minimum height	450 mm
Maximum height	880 mm
Dimensions folded	105x485x65 mm
Materials	Al, Nylon
Weight	560 g

<sup>1</sup> Heights measured from the stretcher frame with opened IV holders

Class I MD according to dir. 93/42/CEE

TRACK 4 IV HOLDER FOR  
506/8409/CROSS-CHAIR  
IF01043B



CND Classification	V80
Registration number	1229046

Rev.0 (10/02/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)