

## TRACK 4

### IV pole for stretcher



Track 4 is an IV pole for Spencer stretcher, intended to maintain the prevalence necessary in order to guarantee the continuity of the infusion operations also during the transport with the stretcher.

#### Specific characteristics

- Simple to install
- Rugged and lightweight design
- It guarantees compact size through the ability to be folded up.
- The locking rings allow a quick and effective adjustment
- It can support two bags
- Bands with locking straps to prevent the fluctuations of the bags during the transport.

#### Stretcher compatibility

- Series CROSS / CROSS OVER (no CROSS CHAIR and CROSS UP)
- Series CARRERA
- Series ENDURO

#### Technical data <sup>1</sup>

Minimum height	450 mm
Maximum height	880 mm
Dimension folded	105x485x65 mm
Material	Al, Nylon
Weight	560 gr±50gr

<sup>1</sup> Heights measured from the stretcher frame and with IV open holders

MD of I class in according with Dir. 93/42/CEE

TRACK 4-30 TELESCOPIC IV  
POLE FOR STRETCHER  
IF01047C



CND classification	V80
Registration identification	484826

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)