

TRACK 4-E

IV pole for stretcher



Track 4-E is an IV pole for Spencer stretcher, intended to maintain the prevalence necessary in order to guarantee the continuity of the infusion operations also during the transport with the stretcher.

Specific characteristics

Simple to install
 Rugged and lightweight design
 It guarantees compact size through the ability to be folded up.
 The locking rings allow a quick and effective adjustment
 It can support two bags
 Bands with locking straps to prevent the fluctuations of the bags during the transport.

Technical data ¹

| | |
|------------------|-------------|
| Minimum height | 450 mm |
| Maximum height | 880 mm |
| Dimension folded | 485x65 mm |
| Material | Al, Nylon |
| Weight | 560 gr±50gr |

¹ Heights measured from the stretcher frame and with IV open holders

MD of I class in according with Dir. 93/42/CEE

TRACK 4 TELESCOPIC IV POLE
 SIDEBAR FIXATION
 IF01048C



| | |
|-----------------------------|-------|
| CND classification | V80 |
| Registration identification | 94515 |

Rev.1 (09/08/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it