

## TRACK 6

### IV-bottle holder 10G for recessed mounting



Device intended to support and protect the containers of medicinal substances administered to the patient by venous infusion. Allows a significant reduction of the oscillations that may damage the containers of liquids.

#### Specific features

- Easy ceiling installation
- Robust and light structure
- 10g fixing of the bottle
- Space saving device
- Can hold 250 or 500cc bottles
- Extremely easy to use

#### Technical data

Dimensions closed*	280x150x40 mm
Dimensions opened*	280x150x180 mm
Recessed part	240x110x20 mm
Material	Nylon
Loading capacity**	0,5 kg
Weight	500 ± 20g

\*Dimensions refer to the installed device

\*\* 500g or 0,5L bottle

#### Compliant with

EN 1789

TRACK 6 TWIN I.V.POLE  
HOLDER  
**IF01050B**

Rev.2 (11/11/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)

© Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA S.r.l.