



TRACK 8

Compact 10G IV-bottle holder



Device intended to support and protect the containers of medicinal substances administered to the patient by venous infusion. Allows a significant reduction of the oscillations that may damage the containers of liquids.

Specific Characteristics

- Easy ceiling installation
- Robust and light structure
- 10g fixing of the bags
- Space saving device
- Can hold 100,250 or 500cc bags
- Practical to use thanks to the locking system and straps

Compliant with

EN 1789

Technical data

| | |
|-------------------|-------------------|
| Dimensions closed | 240x68x45 ± 5 mm |
| Dimensions opened | 215x68x100 ± 5 mm |
| Materials | Nylon, Steel |

Class I Med.Device according to Dir. 93/42/CEE

TRACK 8 IV POLE

IF01052B



CND Classification V80

Registration number 270036

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it