



TRACK 8

Compact 10G IV-bottle holder



Device intended to support and protect the containers of medicinal substances administered to the patient by venous infusion. Allows a significant reduction of the oscillations that may damage the containers of liquids.

Specific Characteristics

Easy ceiling installation
Robust and light structure
10g fixing of the bags
Space saving device
Can hold 100,250 or 500cc bags
Practical to use thanks to the locking system and straps

Compliant with

EN 1789

Technical data

Dimensions closed	240x68x45 ± 5 mm
Dimensions opened	215x68x100 ± 5 mm
Materials	Nylon, Steel

TRACK 8 IV POLE
IF01052B

Rev.2 (11/11/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it

© Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA S.r.l.