

TRACK 9

IV-bottle holder 10G for recessed mounting



Device intended to support and protect the containers of medicinal substances administered to the patient by venous infusion. Allows a significant reduction of the oscillations that may damage the containers of liquids.

Specific features

- Easy ceiling installation
- Robust and light structure
- 10g fixing of the bottle
- Space saving device
- Can hold 250 or 500cc bottles
- Extremely easy to use

Technical data

Dimensions closed*	260x130x10 mm
Dimensions opened*	260x130x180 mm
Recessed part	226x95x25 mm
Material	Aluminum, steel, nylon
Loading capacity**	0,5 kg

*Dimensions refer to the installed device

** 500g or 0,5L bottle

Class I Med.Device according to Dir. 93/42/CEE

TRACK 9 ALUMINUM IV POLE
IF01054B



CND Classification	V80
Registration number	270019

Rev.0 (31/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it