

TRACK 9 IV-bottle holder 10G for recessed mounting





Device intended to support and protect the containers of medicinal substances administered to the patient by venous infusion. Allows a significant reduction of the oscillations that may damage the containers of liquids.

	Specific features
	Easy ceiling installation
	Robust and light structure
	10g fixing of the bottle
	Space saving device
	Can hold 250 or 500cc bottles
	Extremely easy to use

Technical data	
Dimensions closed*	260x130x10 mm
Dimensions opened*	260x130x180 mm
Recessed part	226x95x25 mm
Material	Aluminum, steel, nylon
Loading capacity**	0,5 kg

^{*}Dimensions refer to the installed device

TRACK 9 ALUMINUM IV POLE

IF01054B

^{** 500}g or 0,5L bottle