

TRACK 5

Telescopic IV pole for Enduro



Track V is an IV pole for Spencer Enduro stretcher, intended to maintain the prevalence necessary in order to guarantee the continuity of the infusion operations also during the transport with the stretcher.

Specific features

- Simple to install
- Rugged and lightweight design
- It guarantees compact size through the ability to be folded up.
- The locking rings allow a quick and effective adjustment
- It can support two bags up to 1,5 kg
- Bands with locking straps to prevent the fluctuations of the bags during the transport.
- Aluminum support hooks for maximum strength

Technical data ¹

Minimum height	430 ± 10 mm
Maximum height	860 ± 10 mm
Width folded	500 ± 10 mm
Materials	Al, Nylon
Weight	1000 ± 20 g

¹ Heights measured from the stretcher frame and with IV open holders

Class I Med. Dev. in according with Dir. 93/42/CEE

TRACK 5 TELESCOPIC IV POLE
FOR ENDURO
IF01055



Rev.3 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it