

## TRACK 6 TWIN

### IV-bottle holder 10g for recessed mounting



Device intended to support and protect the containers of medicinal substances administered to the patient by venous infusion. Allows a significant reduction of the oscillations that may damage the containers of liquids.

#### Specific features

Easy ceiling installation  
Robust and light structure  
10g fixing of the bags and bottles  
Space saving device  
Can hold 250 or 500cc bags  
Practical to use thanks to the locking system and straps  
Gasket on the side in contact with the roof to limit vibrations

#### Technical data

Dimensions closed*	280x150x40 mm
Dimensions opened*	280x150x180 mm
Recessed part	240x110x20 mm
Materials	Nylon, Steel
Loading capacity**	1 kg (equal to 10kg retention capacity)
Weight	530 ± 20g

\*Dimensions refer to the installed device

\*\* 500g+500g or 1L bag

#### Compliant with

EN 1789

TRACK 6 TWIN I.V.POLE  
HOLDER IN ALUMINUM  
**IF01060B**

Rev.3 (11/11/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)

© Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA S.r.l.