

TRACK 6 TWIN

IV-bottle holder 10 for recessed mounting



Device intended to support and protect the containers of medicinal substances administered to the patient by venous infusion. Allows a significant reduction of the oscillations that may damage the containers of liquids.

Specific features

- Easy ceiling installation
- Robust and light structure
- 10g fixing of the bags and bottles
- Space saving device
- Can hold 250 or 500cc bags
- Practical to use thanks to the locking system and straps
- Gasket on the side in contact with the roof to limit vibrations

Technical data

Dimensions closed*	280x150x40 mm
Dimensions opened*	280x150x180 mm
Recessed part	240x110x20 mm
Materials	Nylon, Steel
Loading capacity**	1 kg
Weight	530 ± 20g

*Dimensions refer to the installed device

** 500g+500g or 1L bag

Compliant with

EN 1789

Class I Med.Device according to Dir. 93/42/CEE

TRACK 6 TWIN I.V.POLE
HOLDER IN ALUMINUM
IF01060B



CND Classification	V80
Registration number	270184

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it