


TRACK 11

Compact IV-bottle holder



Device intended to support and protect the containers of medicinal substances administered to the patient by venous infusion. Allows a significant reduction of the oscillations that may damage the containers of liquids.

Specific Characteristics

- Easy ceiling or recessed installation
- Robust and light structure
- 10g fixing of the bags
- Space saving device
- Can hold 250 or 500cc bags
- Practical to use thanks to the locking system and straps

Compliant with

EN 1789

Technical data

Dimensions closed	210x110x27 mm
Dimensions opened	50x120x180 mm
Materials	Acciaio, Al
Weight	290 g

TRACK 11 - COMPACT IV-
BOTTLE HOLDER IN
ALUMINIUM
IF01061B

Rev.0 (11/11/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it

© Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA S.r.l.