

TRACK 4

IV pole for rail and



Track 4 is an IV pole for rail and Sherpa, intended to maintain the prevalence necessary in order to guarantee the continuity of the infusion operations also during the transport.

Specific characteristics

- Simple to install
- Rugged and lightweight design
- It guarantees compact size through the ability to be folded up.
- The locking rings allow a quick and effective adjustment
- It can support two bags
- Bands with locking straps to prevent the fluctuations of the bags during the transport.

Stretcher compatibility

Rail bar and SHERPA

Technical data ¹

Minimum height	450 mm
Maximum height	880 mm
Dimension folded	105x485x65 mm
Material	Al, Nylon
Weight	560 gr±50gr

¹ Heights measured from rail frame and with IV open holders

MD of I class in according with Dir. 93/42/CEE

TRACK 4 FOR RAIL/SHERPA

IF01064C

CND classification



V80

Registration identification

403910

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it