

TRACK V

Telescopic IV pole



Track V is a telescopic IV holder intended to maintain the prevalence necessary in order to guarantee the continuity of the infusion operations also during the transport with the stretcher. This version is dedicated to Spencer 506 and 8409 stretchers.

Specific features

- Easy to install
- Rugged and lightweight design
- It guarantees compact size through the ability to be folded up.
- The locking rings allow a quick and effective adjustment
- It can support two bags
- Bands with locking straps to prevent the fluctuations of the bags during the transport.

Compatible stretchers

- Spencer 506
- Cross Up 8409
- Cross Chair

Technical data

Minimum height	410 ± 20 mm
Maximum height	850 ± 20 mm
Dimensions folded	510x95x65 ± 20 mm
Materials	Al, Nylon
Loading capacity	1 kg

Class I med.dev. according with Dir. 93/42/CEE

TRACK 5 TELESCOPIC IV POLE
FOR 506/8409
IF01067B



CND Classification	V80
Registration number	958635

Rev.1 (06/12/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it