

MEGA HOT

Thermos insulating bag



Device intended to sustain and to provide a thermal protection to the content.

Specific Characteristics

- Transparent window on the side for monitoring the amount of solution
- A transparent pocket on the side can be used for labels to identify the content
- On the back there is a male strap bend
- The device is internally made of soft synthetic fabric
- On the top, the belt can be used as handle or as an attachment to another holding device
- The bottle is easily stored thanks to the strap closure

Technical data

Dimensions	320x140x70 ± 10mm
Storable bottles	From 500 to 1000cc
Materials	Nylon, Polyester,PVC
Weight	90 g

Class I Med.Device according to Dir. 93/42/CEE

MEGA HOT - THERMO-INSULATING HEATING BAGS

IF02030B



CND Classification Z12030305

Registration number 304172

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it