

NEC LOC

Cervical collar



Bivalve cervical collar designed for use in first aid and emergency situations. It's identified as a solution for the immobilization of the cervical spine in trauma patients.

Specific Characteristics

Quick application due to strap closure
 Access for tracheotomy
 Comfortable soft padding

Technical data

Size	Pediatric 6-12 years
Materials	PE/Plastizote
Weight	170 g

Class I Med.Device according to Dir. 93/42/CEE

NEC LOC CERVICAL COLLAR
 RIGID 6-12 YEARS
 JM10009A



Distributed by Spencer Italia S.r.l.

Rev.1 (31/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it